



Myelin Basic Protein (MBP)

Concentrated Polyclonal Antibody

Control Number: 901-334-030508

ISO
9001:2000
CERTIFIED

Catalog Number: CP 334 AK, BK
Description: 0.1, 0.5 ml, concentrated
Dilution: 1:500-1:750
Diluent: Van Gogh Yellow

Intended Use:
For In Vitro Diagnostic Use

Summary and Explanation:

Myelin basic protein (MBP) is involved in the myelination of nerves within the central nervous system (CNS). Studies using knockout mice have demonstrated severe neurological deficits such as seizures, tremors and an early death when this protein is lacking. The demyelination disease, Multiple Sclerosis, is an example of a human disease that involves the selective loss of MBP by autoimmune mechanisms. Individuals afflicted with MS develop a compromised neurological condition. Biocare's antibody to MBP can be used to identify alterations or changes in the distribution of this protein in tissue sections.

Principle of Procedure:

Antigen detection, in tissues and cells, is a multi-step immunohistochemical process. The initial step binds the primary antibody to its specific epitope. After labeling the antigen with a primary antibody, a universal, affinity-purified, secondary antibody is added to bind to the primary antibody. An enzyme label is then added to bind to the secondary antibody; this detection of the bound antibody is evidenced by a colorimetric reaction.

Source: Rabbit polyclonal

Species Reactivity: Human

Clone: N/A

Isotype: N/A

Epitope/Antigen: Myelin Basic Protein (MBP)

Cellular Localization: Cytoplasmic, membrane

Positive Control: Brain

Normal Tissue: Brain

Abnormal Tissue: Neurodegenerative tissues

Total Protein Concentration: ~10 mg/ml

Known Applications:

Immunohistochemistry (formalin-fixed paraffin-embedded tissues)

Supplied As: Buffer with protein carrier and preservative.

Van Gogh Yellow Diluent (PD902)

Storage and Stability:

Store at 2°C to 8°C. Do not use after expiration date printed on vial. If reagents are stored under conditions other than those specified in the package insert, they must be verified by the user. Diluted reagents should be used promptly; any remaining reagent should be stored at 2°C to 8°C.

Protocol Recommendations

Peroxide Block:

If using an HRP system, block for 5 minutes with BIOCARE's PEROXIDAZED 1.

Pretreatment Solution (recommended): Diva

Pretreatment Protocol:

Heat Retrieval Method:

Retrieve sections at 80°C for 1 hour using BIOCARE's Decloaking Chamber, followed by a wash in distilled water. Alternatively, steam tissue sections for 45-60 minutes. Allow solution to cool for 20 minutes then wash in distilled water. A significant increase in titer is achieved using Biocare's Diva retrieval solution.

Protein Block:

Incubate for 10-15 minutes at RT with BIOCARE's Background Sniper.

Primary Antibody: Incubate for 60 minutes at RT.

Probe: N/A

Polymer: Incubate for 30 minutes at RT with a Polymer.

Chromogen:

Incubate for 5 minutes at RT when using BIOCARE's DAB. - OR - Incubate for 10 minutes at RT when using BIOCARE's Vulcan Fast Red.

Technical Note:

This antibody has been standardized with BIOCARE's MACH 2 detection system. It can also be used on an automated staining system and with other BIOCARE polymer detection kits. Use TBS buffer for washing steps.

Performance Characteristics:

The optimum antibody dilution and protocols for a specific application can vary. These include, but are not limited to: fixation, heat-retrieval method, incubation times, tissue section thickness and detection kit used. Due to the superior sensitivity of these unique reagents, the recommended incubation times and titers listed are not applicable to other detection systems, as results may vary. The data sheet recommendations and protocols are based on exclusive use of BIOCARE products. Ultimately, it is the responsibility of the investigator to determine optimal conditions. These products are tools that can be used for interpretation of morphological findings in conjunction with other diagnostic tests and pertinent clinical data by a qualified pathologist.

Quality Control:

Refer to NCCLS Quality Assurance for Immunocytochemistry approved guidelines, December 1999 MM4-A Vol.19 No.26 for more information about Tissue Controls.

Precautions:

This antibody contains less than 0.1% sodium azide. Concentrations less than 0.1% are not reportable hazardous materials according to U.S. 29 CFR 1910.1200, OSHA Hazard communication and EC Directive 91/155/EC.

Sodium azide (NaN₃) used as a preservative is toxic if ingested. Sodium azide may react with lead and copper plumbing to form highly explosive metal azides. Upon disposal, flush with large volumes of water to prevent azide build-up in plumbing. (Center for Disease Control, 1976, National Institute of Occupational Safety and Health, 1976)

Specimens, before and after fixation and all materials exposed to them, should be handled as if capable of transmitting infection and disposed of with proper precautions. Never pipette reagents by mouth and avoid contacting the skin and mucous membranes with reagents and specimens. If reagents or specimens come in contact with sensitive areas, wash with copious amounts of water.

Microbial contamination of reagents may result in an increase in nonspecific staining. Incubation times or temperatures other than those specified may give erroneous results. The user must validate any such change. The MSDS is available upon request.

Troubleshooting:

Follow the antibody specific protocol recommendations according to data sheet provided. If atypical results occur, contact BIOCARE's Technical Support at 1-800-542-2002.

Limitations and Warranty:

There are no warranties, expressed or implied, which extend beyond this description. BIOCARE is not liable for property damage, personal injury, or economic loss caused by this product.

References:

1. Reindl M et al. Antibodies against the myelin oligodendrocyte glycoprotein and the myelin basic protein in multiple sclerosis and other neurological diseases: a comparative study. Brain. 1999 Nov; 122 (Pt 11):2047-56.
2. Johnson MD, Glick AD, Davis BW. Immunohistochemical evaluation of Leu-7, myelin basic-protein, S1000-protein, glial-fibrillary acidic-protein, and LN3 immunoreactivity in nerve sheath tumors and sarcomas. Arch Pathol Lab Med. 1988 Feb;112(2):155-60.





Myelin Basic Protein (MBP)

Concentrated Polyclonal Antibody

Control Number: 901-334-030508

ISO
9001:2000
CERTIFIED**References cont'd:**

3. Natalia A. Ponomarenko et al. Autoantibodies to myelin basic protein catalyze site-specific degradation of their antigen. PNAS 2006 Jan;103 (2):281-86.
4. Center for Disease Control Manual. Guide: Safety Management, NO. CDC-22, Atlanta, GA. April 30, 1976 "Decontamination of Laboratory Sink Drains to Remove Azide Salts."
5. National Committee for Clinical Laboratory Standards(NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. Villanova, PA 1991;7(9). Order code M29-P.

